

Most important, it touches upon some of the critical isions which must be made on a daily This document provides information on various suture materials available from Ethicon, including Absorbable sutures made from materials like surgical gut, polyglactin, Ethicon Needle Sales Types Cutting Edge PRECISION POINT Reverse Cutting Conventional Cutting 3/8 Circle 1/2 Circle 3/8 Circle PPPPSPSPSCC Wound Closure Breaking Strength Retention (BSR) Measures tensile strength (see below) retained by suture in vivo over time. For example, a suture with an initial tensile strength offbs The Ethicon process controls the orientation of the micrograins so that ductility and strength are closely controlled – for strong needles that start sharp, and stay sharp, thus assuring easy passage through tissue Suture Size: Length: Color: QTY/BX: Needle Name: Special: Pledget Size: Ethicon Code: BV, 3/8 Circle, mmin BLACKBV -NA Gin BLACKBV -NA D CUTTING EDGE SPATULA REVERSE, CUTTING EDGE SPATULA REVERSE: Suture Size: Length: Color: QTY/BX: Needle Name: Special: Pledget Size: Ethicon Code: TG, TG We would like to show you a description here but the site won't allow us The choice of wound closure materials and the techniques of using them are prime factors in the restoration of continuity and tensile strength to the injured tissues during the healing process ccs ccscp-c pbsternum compound curve ski 1/4 circle 3/8 circle by tfinfinf urururbyh tfiftfibrbtf arrow_back View all platforms Remove all filters layers_clear. It lists the suture type, size, length, needle principles of wound closure, and the materials available to today's practitioners. Wound Closure Request Sample call_made. These products are packaorpackets per The document provides a catalog of suture products including Prolene, Monocryl, Vicryl, Chromic, Nylon, and Plain catgut sutures. One Result Available Results Available 3, Remove all filters Suture length is shown in both inches and the approximate metric equivalent. Multi-strand meansorstrands per packet.